

Lester Proximal Tibial TKA Retractor

Designed by D. Kevin Lester, MD

Helps expose the cut surface of the tibia to allow sizing, preparation and cleansing during TKA

Also helps protect the posterior knee soft tissue structures from injury.



PRODUCT NO:
4699
Overall Length: 12" (30,5 cm)
Depth from Bend: 5" (12,7 cm)
Blade Width: 1.5" (3,8 cm)



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

INNOMED



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362

© 2023 Innomed, Inc.