

Lighted Fukuda-type Retractors

Designed by Evan Flatow, MD & Louis Bigliani, MD

Used to retract the humeral shaft posteriorly and helping to expose the entire glenoid surface, the lighting attachment helps provide enhanced visual exposure

Can be attached to a fiber optic light cable with ACMI (female) connector.

Lighted retractors come with one (1) Disposable LED Light Source (#8010-01).

PRODUCT NO'S:	
1930-L-01 [Lighted Narrow]	
Blade Width: 32 mm	
Opening: 25 x 40 mm	
Overall Length: 8.75" (22,2 cm)	
1940-L-01 [Lighted Wide]	
Blade Width: 38 mm	
Opening: 32 x 40 mm	
Overall Length: 8.75" (22,2 cm)	



Light Source Cable Adapters

PRODUCT NO'S:	
8009-S [ACMI to Storz Adapter]	
8009-W [ACMI to Wolf Adapter]	



Disposable LED Light Source



PRODUCT NO'S:	
PACKAGE OF 1:	
8010-01 [Disposable LED Light Source]	
Overall Length: 2.5" (6,4 cm)	
Diameter: 1" (2,54 cm)	
PACKAGE OF 10:	
8010-10 [Disposable LED Light Source]	

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



INNOMED



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362