

Lighted Hip Retractors

Lighting attachment for enhanced visual exposure

Excellent for use in acetabular exposure and total hip replacements. Especially useful for anterior approach.

Can be attached to a fiber optic light cable with ACMI (female) connector. Retractors can be steam sterilized.



*Jana Lighted
Cobra Retractors*

Designed by Ajoy K. Jana, MD

Lighted Cobra Retractors

*Lighted Single Prong Double
Bent Hohmann Acetabular Retractor - Long*

PRODUCT NO'S:
6119-L [Jana Lighted Cobra] Overall Length: 14.75" (37,5 cm) Blade at Widest: 33 mm
6120-L [Lighted Narrow Cobra] Overall Length: 11.75" (29,8 cm) Handle Length: 6.5" (16,5 cm) Blade Width: 19 mm
6130-L [Lighted Standard Cobra] Overall Length: 12" (30,5 cm) Handle Length: 7" (17,8 cm) Blade Width at Widest: 33 mm
6135-L [Lighted Deep Cobra] Overall Length: 16" (40,6 cm) Handle Length: 7" (17,8 cm) Blade Width at Widest: 33 mm
6210-02L [Lighted Single Prong Double Bent Hohmann - Long] Overall Length: 12.5" (31,8 cm) Blade and Tip Length: 3" (76 mm) Blade Width: 15 mm



Light Source Cable Adapters

PRODUCT NO'S:
8009-S [ACMI to Storz Adapter]
8009-W [ACMI to Wolf Adapter]



www.innomed.net

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

INNOMED



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362

© 2023 Innomed, Inc.