

Lighted Hip Retractors


Lighting attachment for enhanced visual exposure

Excellent for use in acetabular exposure and total hip replacements. Especially useful for anterior approach.

Can be attached to a fiber optic light cable with ACMI (female) connector. Retractors can be steam sterilized.

Jana Lighted Cobra Retractors

Designed by Ajoy K. Jana, MD



Lighted Cobra Retractors

Lighted Single Prong Double Bent Hohmann Acetabular Retractor - Long

PRODUCT NO'S:
6119-L [Jana Lighted Cobra] Overall Length: 14.75" (37,5 cm) Blade at Widest: 33 mm
6120-L [Lighted Narrow Cobra] Overall Length: 11.75" (29,8 cm) Handle Length: 6.5" (16,5 cm) Blade Width: 19 mm
6130-L [Lighted Standard Cobra] Overall Length: 12" (30,5 cm) Handle Length: 7" (17,8 cm) Blade Width at Widest: 33 mm
6135-L [Lighted Deep Cobra] Overall Length: 14.5" (36,9 cm) Handle Length: 7" (17,8 cm) Blade Width at Widest: 33 mm
6210-02L [Lighted Single Prong Double Bent Hohmann - Long] Overall Length: 12.5" (31,8 cm) Blade and Tip Length: 3" (76 mm) Blade Width: 15 mm



Light Source Cable Adapters

PRODUCT NO'S:
8009-S [ACMI to Storz Adapter]
8009-W [ACMI to Wolf Adapter]



www.innomed.net

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

INNOMED



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362

© 2023 Innomed, Inc.