

Lighted Yankaur Suction Device

Designed by Adolph V. Lombardi Jr., MD

Designed to help provide effective suction with the addition of a light source for enhanced visualization



- ▶ Device comes with one (1) Disposable LED Light Source (#8010-01)
- ▶ Can also be attached to a fiber optic light cable with ACMI (female) connector
- ▶ The handle is made of Delrin
- ▶ Entire device is steam sterilizable

Disposable LED Light Source



PRODUCT NO'S:

PACKAGE OF 1:

8010-01 [Disposable LED Light Source]
Overall Length: 2.5" (6,4 cm)
Diameter: 1" (2,54 cm)

PACKAGE OF 10:

8010-10 [Disposable LED Light Source]

PRODUCT NO'S:


8016-L-01 [Kit – Device with Disposable LED Light Source]

Overall Length: 11.75" (29,8 cm)
Handle Length: 3.93" (10 cm)
Handle Width: .86" (2,2 cm)
Suction Tube Diameter: .25" (6,35 mm)

Also available individually:

8016-L [Suction Device Only]

Overall Length: 11.75" (29,8 cm)
Handle Length: 3.93" (10 cm)
Handle Width: .86" (2,2 cm)
Suction Tube Diameter: .25" (6,35 mm)



ISO 13485:2016

Scan to
Launch Our
Website



FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

INNOMED

1.800.548.2362

© 2019 Innomed, Inc.