



Lighted Yankaur Suction Device

Designed by Adolph V. Lombardi Jr., MD

Designed to help provide effective suction with the addition of a light source for enhanced visualization



Can be attached to a fiber optic light cable with ACMI (female) connector.
The handle is made of Delrin. Entire device is steam sterilizable.

PRODUCT NO

8016-L

Overall Length: 11.75" (29.8 cm)

Handle Length: 3.93" (10 cm)

Handle Width: .86" (2.2 cm)

Suction Tube Diameter: .25" (6.35 mm)



USA MADE

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

INNOMED

1.800.548.2362