

Standard Handle

Lombardi Femoral Tibial Spreader

Designed by Adolph V. Lombardi Jr., MD

Thin pads help to separate the femur and tibia during total knee procedures



The calibrated handle of the spreader helps to accurately gauge the gap, and makes it possible for two spreaders to be used to assist in balancing ligaments.



Large



Small



Large with Locking Mechanism

Locking ratchet mechanism helps prevent accidental release, and provides for controlled adjustment and easy release.



Small with Locking Mechanism

Designed with the grip closer together for easier gripping and to help reduce hand fatigue

New!



Small with Small Grip

PRODUCT NO'S:

1875 [Large – Grooved Pads]

Overall Length: 9.25" (23,5 cm)

Pads: 22 mm x 13 mm

Opens to: 50 mm

1875-D [Large – Diamond Cut Pads]

Overall Length: 9.25" (23,5 cm)

Pads: 22 mm x 13 mm

Opens to: 50 mm

1876 [Small – Grooved Pads]

Overall Length: 7" (17,8 cm)

Pads: 22 mm x 13 mm


Opens to: 35 mm

1876-D [Small – Diamond Cut Pads]

Overall Length: 7" (17,8 cm)

Pads: 22 mm x 13 mm

Opens to: 35 mm



Grooved Pads Diamond Cut Pads

Locking Mechanism

PRODUCT NO'S:

1875-LR [Large – Grooved Pads with Locking Mechanism]

Overall Length: 9.625" (24,4 cm)

Pads: 22 mm x 13 mm


Opens to: 45 mm

1876-LR [Small – Grooved Pads with Locking Mechanism]

Overall Length: 7" (17,8 cm)

Pads: 22 mm x 13 mm

Opens to: 35 mm



Small Grip Handle

PRODUCT NO:

1876-SG [Small with Small Grip and Grooved Pads]

Overall Length: 7" (17,8 cm)

Pads: 22 mm x 13 mm

Opens to: 35 mm



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

INNOMED 



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362

© 2021 Innomed, Inc.