



# Lombardi Femoral Tibial Spreader with Easy Release Locking Mechanism

Spreader designed by Adolph V. Lombardi Jr., MD.  
Locking mechanism designed by Munish C. Gupta, MD

*Thin pads help to separate the femur and tibia during total knee procedures*

Locking ratchet mechanism helps prevent accidental release, and provides for controlled adjustment and easy release.



**PRODUCT NO'S:**  
1875-LR  
Overall Length: 9.625" (24,4 cm)  
Pads: 22 mm x 13 mm  
Opens to 45 mm



ISO 13485:2016

**FREE TRIAL ON MOST INSTRUMENTS**

Scan to  
Launch Our  
Website



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
Fax +41 41 740 67 71

**INNOMED**

**1.800.548.2362**