




Lombardi Gap Balancing Femoral Tibial Spreader with Easy Release Locking Mechanism

Spreader designed by Adolph V. Lombardi, Jr., MD.
Locking mechanism designed by Munish C. Gupta, MD

Designed to help separate the femur and tibia during total knee procedures, with the pads being parallel when measured at 20mm of separation

Locking ratchet mechanism helps prevent accidental release, and provides for controlled adjustment and easy release.



PRODUCT NO:
1878-LR
Overall Length: 9.625" (24,4 cm)
Pads: 22 mm x 13 mm
Opens to 45 mm



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

INNOMED 

1.800.548.2362