



Lombardi Hip Cup Liner/Shell Extractor

Designed by Adolph V. Lombardi, MD

Used for removal of a total hip cup or liner

Expandable flanges are designed to bite into the polyethylene of a total hip cup. When the flanges have been expanded, a slap hammer is screwed into the extractor for removal. The extractor can also be used for removal of a metal hip cup shell if the shell has a groove around the rim for the flanges to lock into. Also very helpful for cemented cup extraction. Set includes standard slap hammer #3925.



PRODUCT NO:
3638-00 [Set]
<i>Also Available Individually</i>
3638-01 [Remover Only]
<i>Overall Length: 9.5" (24,1 cm)</i>
3925 [Standard Slap Hammer]
<i>3/8"-16 Thread Gauge</i>

**PROUDLY
MADE
IN THE
USA**



ISO 9001:2008 • ISO 13485:2003

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



INNOMED



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362

© 2015 Innomed, Inc.