



Lombardi Leg Positioner

Designed by Adolph V. Lombardi Jr., MD

Designed to hold the leg during total knee surgery, the unrestricted design helps allow for manipulation of the leg

- ▶ The footpiece consists of an open topped boot, separated with four spacers and attached to a boot support plate.
- ▶ The support plate is bent beyond the rear of the boot to help provide support in different positions of flexion/extension.
- ▶ The spacers allow for room to clean between the boot and support plate.



Leg Positioner Sterile Protective Pad & Wrap

Disposable, latex-free sterile foam pad and cohesive wrap helps protect patient from pressure sores, abrasions and possible neurological impairment while securing foot into the boot.

PRODUCT NO'S:
2629-00 [Case of 10 Sets]
2629-L [1 Set]



Each set includes 1 Pad and 1 Wrap



PRODUCT NO'S:
2622 [Lombardi Leg Positioner Kit]
Kit Includes:
2622-FP [Footpiece Only] Dimensions: 13.3" x 12.8" x 7.5" (33,8 x 32,5 x 19 cm)
2629-L [Sterile Pad & Wrap - 1 Set] Two included in kit; one with this product number
Replacement Parts:
2629-00 [Sterile Pad & Wrap - Case of 10 Sets]



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

INNOMED



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362

© 2021 Innomed, Inc.