



Prototype Shown

Lubahn Carpal Corkscrews

Designed by John D. Lubahn, MD

Designed to help with removal of carpal and/or tarsal bones

- ▶ Aids trapezium removal during basal joint arthroplasty when the bone is being removed as a unit
- ▶ Can also be used to facilitate a proximal row carpectomy as it fits the scaphoid, lunate, and triquetrum
- ▶ May additionally be used to remove the pisiform in cases of arthritis of the piso-triquetral joint



PRODUCT NO'S:	
1191	[Standard] Overall Length: 2.25" (5,7 cm)
1191-01	[Extended] Overall Length: 6.5" (16,5 cm)



Prototype Shown

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



INNOMED 



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362