




Minimally Invasive Knee Retractors

Helps provide excellent visibility and ligament protection during Total and Unicondylar Knee Replacement Surgery


PRODUCT NO'S:	
1	S3035 [Small Hohmann Retractor] Overall Length: 7.5" (19,1 cm) Blade Width: 25 mm
2	S3037 [Condylar Retractor] Overall Length: 7.5" (19,1 cm) Blade Width: 12 mm
3	S3038 [Superior Retractor] Overall Length: 8.75" (22,2 cm) Blade Width: 31 mm
4	S3039 [Patellar Retractor] Overall Length: 10.25" (26 cm) Blade Width: 45 mm



1 MI Small Hohmann Retractor



3 MI Superior Retractor



2 MI Condylar Retractor



4 MI Patellar Retractor

Retractor handles are designed for optional use with modular weights to help hold the retractors in place.



Modular Weights

PRODUCT NO'S:	
3430-01	1.5 lbs. (.68 kg)
3430-02	2.0 lbs. (.91 kg)
3430-03	2.5 lbs. (1.13 kg) with attaching hook

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

www.innomed.net



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362

© 2025 Innomed, Inc.