



# Macko Square Tipped Rongeur

Designed by Victor W. Macko, MD

*Square tipped rongeur with an ergonomic grip, double action mechanism, long reach, and low profile – for use in total ankle, hip, knee, and spine surgery*

In total ankle surgery, the aggressive, low profile jaws aid in the removal of tibial bone in spite of limited space. The square ended tip helps produce a flat, finished surface following anterior talar facet reaming.

When used for morcelizing bone graft, the shallow, wide jaw helps avoid impaction.



Three Jaw Sizes Available

PRODUCT NO'S:
1778-01 Jaw Bite: 5 x 18 mm Overall Length: 10" (25,4 cm)
1778-02 Jaw Bite: 7 x 18 mm Overall Length: 10" (25,4 cm)
1778-03 Jaw Bite: 10 x 18 mm Overall Length: 10" (25,4 cm)



[www.innomed.net](http://www.innomed.net)

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

# INNOMED



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
info@innomed-europe.com

1.800.548.2362