

# Malleable Bone Tamp – Extra Small

Modified by Serge Kaska, MD & Amal Das, MD


*Designed to help impact bone into acetabular cup holes*

**PRODUCT NO:**

5296-02 [Extra Small]  
Overall Length: 11.4" (29cm)  
Shaft Length: 5.9" (15cm)  
Impactor Diameter: 6.5mm

**PROUDLY  
MADE  
IN THE  
USA**

*Malleable shaft can be contoured for different angles*



ISO 9001:2008 • ISO 13485:2003

**FREE TRIAL ON MOST INSTRUMENTS**

Scan to  
Launch Our  
Website



**INNOMED**



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

**Innomed-Europe** Tel. +41 41 740 67 74  
Fax +41 41 740 67 71

**1.800.548.2362**