


Malleable Bone Tamps

Modified by Serge Kaska, MD

The large tamp is designed to help elevate a depressed tibial plateau fracture, while the small tamp can help elevate a depressed tibial plafond and smaller tibial plateau fractures

Malleable shaft can be contoured for different angles.

Bend the shaft to get just the angle you need





PRODUCT NO'S:

5296 [Large]
Overall Length: 14" (35,6cm)
Shaft Length: 9.5" (24,1cm)
Impactor Diameter: 12.5mm

5296-01 [Small]
Overall Length: 9.5" (24,1cm)
Shaft Length: 6" (15,2cm)
Impactor Diameter: 10mm

PROUDLY
MADE
IN THE
USA



ISO 9001:2008 • ISO 13485:2003

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Scan to
Launch Our
Website



INNOMED

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362