


Manzary Proximal Tibial Stabilizing Knee Retractor

Design modification by Mojib Manzary, MD, FRCS of original design by D. Kevin Lester, MD, and Christopher M. Meckel, MD

Designed to help subluxe the tibia anteriorly in posterior stabilizing total knee replacement, helping to expose the proximal surface for preparation



PRODUCT NO:
4531
Overall Length: 12.75" (32.4 cm)
Blade Width: 1.5" (3.8 cm)
Block Dimensions: .8" x 1.5" x .5"
(2 x 3.8 x 1.25 cm)



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



INNOMED



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362