


# Mazzara Rongeur with Small Pistol Grip Handle

Designed by James T. Mazzara, MD

*Designed for bone and soft tissue removal in small joint surgery, the small pistol grip handle lessens hand fatigue and slippage, and allows for better visualization*



PRODUCT NO'S:
1765-04 Jaw Bite: 2 x 10 mm Overall Length: 9" (22,9 cm)
1765-05 Jaw Bite: 4 x 10 mm Overall Length: 9" (22,9 cm)



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to  
Launch Our  
Website



**INNOMED** 



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
Fax +41 41 740 67 71

1.800.548.2362