

McMaster Abductor Retractor

Designed by William D. McMaster, MD

Designed to help with proximal femur exposure helping to protect the abductors – gluteus medius and minimus – during posterior approach THA

The ergonomic design allows any application where soft tissue retraction is needed for exposure similar to a doyen or hibbs-type retractor.



PRODUCT NO:
6385
Overall Length: 8.75" (22.2 cm)
Handle Length: 4.625" (11.7 cm)
Blade Depth: 4.125" (10.8 cm)
Blade Width: 1.5" (3.8 cm)
Blade Curve Diameter: 1.75" (4.4 cm)



Profile View

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe

Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362

INNOMED