

Medial Acetabular Retractors

Designed for acetabular exposure during total hip surgery



PRODUCT NO'S:
C1001 [Left] Overall Length: 15" (38,1 cm) Handle Length: 6" (15,2 cm) Blade Width: 18 mm
C1002 [Right] Overall Length: 15" (38,1 cm) Handle Length: 6" (15,2 cm) Blade Width: 18 mm



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



INNOMED



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362