



Meniscal Clamp

Redesigned clamp is curved for easier use, visualization, and tissue holding

PRODUCT NO:


1883

Overall Length: 7" (17,8cm)

Teeth Length: .082" (2mm)

Jaw Length: 1.5" (3,8cm)

**PROUDLY
MADE
IN THE
USA**



ISO 9001:2008 • ISO 13485:2003

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



INNOMED 



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362