


# Merchant Surgical Goniometer

Modified by Al Merchant, MD

*Designed to help assess frontal plane limb alignment or measure the Q angle*

The extended length can reach from the center of the knee to the femoral head or the anterior superior iliac spine. The collapsible stainless steel device is autoclavable.



Collapsible and steam sterilizable.

**PRODUCT NO:**

2029  
Overall Length:  
41" Fully Extended (104,2 cm)  
22,5" Folded in Half (57,2 cm)  
12" Fully Collapsed (30,5 cm)



ISO 9001:2008 • ISO 13485:2003

**FREE TRIAL ON MOST INSTRUMENTS**

Scan to  
Launch Our  
Website



**INNOMED** <sup>®</sup>



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
Fax +41 41 740 67 71

**1.800.548.2362**