


# Meyerding Type Retractors with Ergonomic Handle

*Designed for general use soft tissue retraction, the ergonomic handle allows for a better grip and less fatigue*

Non-glare finish featured on the metal retractor parts.



PRODUCT NO:
6241 [50 x 16 mm] Overall Length: 8.875" (22,5 cm) Blade Width: 16 mm Blade Depth: 50 mm
6242 [75 x 15 mm] Overall Length: 9" (22,9 cm) Blade Width: 15 mm Blade Depth: 75 mm
6243 [75 x 25 mm] Overall Length: 9" (22,9 cm) Blade Width: 25 mm Blade Depth: 75 mm



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to  
Launch Our  
Website



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
Fax +41 41 740 67 71

**INNOMED**

1.800.548.2362