


# Micro Curettes

Four cup sizes, straight or 45° angled-end shaft



Ultra hard titanium nitride coating helps to extend curette life by increasing surface hardness, prolonging sharpness, and resisting chemicals and corrosion.

PRODUCT NO'S:	
Straight Micro Curettes	Angled Micro Curettes
Overall Length: 9.75" (24,8 cm) Shaft Length: 4.5" (11,4 cm)	Overall Length: 9.75" (24,8 cm) Shaft Length: 4.5" (11,4 cm)
4242 <b>Cup Size 2</b>	4242-01 <b>Cup Size 2</b>
4240 <b>Cup Size 1</b>	4240-01 <b>Cup Size 1</b>
4244 <b>Cup Size 4/0</b>	4244-01 <b>Cup Size 4/0</b>
4246 <b>Cup Size 6/0</b>	4246-01 <b>Cup Size 6/0</b>



[www.innomed.net](http://www.innomed.net)

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

**INNOMED** 



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
info@innomed-europe.com

1.800.548.2362

© 2023 Innomed, Inc.