



# Modified Lambotte Cup Removal Osteotomes

*Designed with different hemisphere of curves to match cups of different sizes*

Four osteotomes with different hemispherical radii allow the osteotomes to fit next to the outer surface of different size acetabular hip cups. The osteotomes have a handle for better control and a hammering platform.



**PRODUCT NO.'S:**

5240-44

Blade Radius: 44mm  
Overall Length: 12.75"  
Handle Length: 4.75"

5240-48

Blade Radius: 48mm  
Overall Length: 12.75"  
Handle Length: 4.75"

5240-52

Blade Radius: 52mm  
Overall Length: 12.75"  
Handle Length: 4.75"

5240-56

Blade Radius: 56mm  
Overall Length: 12.75"  
Handle Length: 4.75"

**FREE TRIAL ON MOST INSTRUMENTS**



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax



# INNOMED

**1.800.548.2362**