


Modified Anterior Hip Retractor

Designed to provide exposure of the proximal femur

PRODUCT NO:

6421 [Narrow Tip]
Overall Length: 15.75" (40 cm)
Blade Width: 1.15" (3 cm)

6422 [Wide Tip]
Overall Length: 15.75" (40 cm)
Blade Width: 1.15" (3 cm)



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

INNOMED

1.800.548.2362