


Modified Anterolateral Retractor Set

Used for anterior MIS hip surgery



PRODUCT NO'S:

6161-00 [Modified Anterolateral Retractor Set]
Set includes: (2) 6162, (1) 6163, (1) 6164

Also available individually:

6162 [Modified Deep Hohmann Retractor – A]
(2) included in set, (1) only with this product number
Overall Length: 14.5" (36,9 cm)
Blade Width: 25 mm

6163 [Modified Small Hohmann Retractor – B]
Overall Length: 8.5" (21,6 cm)
Blade Width: 18 mm

6164 [Modified Mueller Retractor – C]
Overall Length: 15.25" (38,8 cm)
Blade Width: 25 mm



ISO 9001:2008 • ISO 13485:2003

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

INNOMED

1.800.548.2362