



# Modified Kerrison Punch

*Designed with an ergonomic pistol grip handle and a unique toggle mechanism that allows for a 10:1 jaw force ratio*

- ▶ 40° upward jaw design
- ▶ Disassembles for cleaning
- ▶ Internal bone ejector
- ▶ Four sizes: 2, 3, 4, and 5 mm



PRODUCT NO'S:	
3626 [2 mm Punch]	Overall Length: 12.5" (31,8 cm) Shaft Length: 9" (22,9 cm)
3627 [3 mm Punch]	Overall Length: 12.5" (31,8 cm) Shaft Length: 9" (22,9 cm)
3628 [4 mm Punch]	Overall Length: 12.5" (31,8 cm) Shaft Length: 9" (22,9 cm)
3629 [5 mm Punch]	Overall Length: 12.5" (31,8 cm) Shaft Length: 9" (22,9 cm)



Disassembles for cleaning

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to  
Launch Our  
Website



**INNOMED**



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
Fax +41 41 740 67 71

1.800.548.2362