

Modified Lambotte Osteotomes

Designed with a striking platform, plus a cross-bar hole to help control rotational stability and assist with removal

- ▶ Six (6) sizes available, from 1/4" to 1-1/2" in 1/4" increments.
- ▶ Cross-bar and case included in complete set.
- ▶ Two smallest sizes have an 1/8" hole in which an 1/8" pin can be used as a cross bar (not included).



PRODUCT NO'S:

5350-00 [Set w/Case]

Also Available Individually:

5350-25* [1/4"]

Overall Length: 9" (22,9 cm)
Osteotome Width: .25" (6,4 mm)

5350-50* [1/2"]

Overall Length: 9" (22,9 cm)
Osteotome Width: .5" (12,7 mm)

5350-75 [3/4"]

Overall Length: 9" (22,9 cm)
Osteotome Width: .75" (19 mm)

5350-100 [1"]

Overall Length: 9" (22,9 cm)
Osteotome Width: 1" (25,4 mm)



5350-125 [1-1/4"]

Overall Length: 9" (22,9 cm)
Osteotome Width: 1.25" (31,8 mm)

5350-150 [1-1/2"]

Overall Length: 9" (22,9 cm)
Osteotome Width: 1.5" (38,1 mm)

5350-CASE [Case]

Dimensions: 12.25" x 11.25" x 1"
(31,1 x 28,6 x 2,5 cm)

5350-CB [Cross Bar]

Overall Length: 4.375" (11,1 cm)

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



INNOMED



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362