



Modular Head Holder

Designed by Byron E. Dunaway, MD & Wayne Goldstein, MD

Designed to hold 22mm to 36mm heads for ease of insertion in minimally invasive THR



Head holding ends are plastic coated to help eliminate any damage to the implant. Available in two lengths. Steam and gas sterilizable.

PRODUCT NO'S:
8290-01 [7"] Overall Length: 7" (17,8 cm)
8290-02 [9"] Overall Length: 9" (22,9 cm)



ISO 9001:2008 • ISO 13485:2003

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

INNOMED

1.800.548.2362