



Monaco Small Space Retractor

Designed modified by Spencer Monaco, DPM, FACFAS

Designed to retract adipose tissue and surrounding soft tissue structures through a small incision for open plantar fasciotomies, neuroma excisions and the lateral release during bunion surgery
Also useful for various hand surgeries such as open carpal tunnel surgery.



PRODUCT NO:
1887-01
Overall Length: 4.25" (10,8 cm)
Blade Depth: 18 mm
Blade Width: 12 mm
Blade Lip: 3.5 mm

MADE EXCLUSIVELY
FOR INNOMED IN
GERMANY

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362

INNOMED