

Mongold Capsule Knife

Designed by Evie Mongold, MD

*Designed to reach behind
the femoral head to release
the capsule ligament*


PRODUCT NO:

4115

Overall Length: 7.75" (19,7 cm)

Blade Diameter: 2" (5,1 cm)

Blade Width: .5" (1,3 cm)



www.innomed.net

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

INNOMED



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362

© 2023 Innomed, Inc.