

Mueller-type Femoral Neck Elevator

Designed to elevate the proximal femur

PRODUCT NO:

3415

Overall Length: 13.5"


Handle Length: 6.5"

Blade Width at Widest: 25mm

PROUDLY

MADE
IN THE

USA



ISO 9001:2008 • ISO 13485:2003

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

INNOMED

1.800.548.2362