

Namba Hip Slide

Designed by Robert S. Namba, MD


Gently slides the femoral head into the cup without scratching – essential for ceramic heads

Helps reduce a femoral head trial and implant into the acetabulum during total hip surgery, and helps avoid injury to the knee from excessive traction during hip reduction. Manufactured of delrin to help eliminate damage to the implant. Can be steam or gas sterilized and is radiolucent.



PRODUCT NO'S:	
Overall Length: 12"	
6890	For 22-40mm heads
6891	For 40-48mm heads
6892	For 50-60mm heads

Facilitates MIS hip replacement procedures



FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax



INNOMED 

1.800.548.2362