

1-800-548-2362

International (912) 236-0000



www.innomed.net

(912) 236-7766 Fax

©2002 Innomed, Inc., All Rights Reserved.

**Innomed**  
Instrument Division  
103 Estus Drive, Savannah, GA 31404




FREE TRIAL ON MOST INSTRUMENTS  
CALL FOR MORE INFORMATION

## Nicholson Headrest



Designed to provide excellent exposure to the shoulder, the headrest can be used with standard OR tables (with no modifications to the table). The headrest provides patient support and helps position the patient for all types of shoulder surgery—arthroscopic and open—in the beachchair position. It can be quickly placed and adjusted. A velcro head strap is included for optional use.

Designed by  
Greg Nicholson, M.D.

Product No: 2450



*Helps position the patient for  
all types of shoulder surgery in  
the beachchair position*



Provides excellent exposure to the shoulder, and patient support.