

# O'Reilly Dual Handle Direct Anterior Retractor

Designed by Michael P. O'Reilly, MD

*Designed for use over the anterior pelvic rim during acetabular exposure in direct anterior THA, the dual handle design allows for comfortable use in both right and left hips, as well as easy exchange of the instrument between assistants*

Can be used in MIS/Direct Anterior, Total Hip Arthroplasty, Posterior/Anterolateral THA, and Hemiarthroplasty.

**PRODUCT NO:**

3011  
Overall Length: 13.25"  
Blade Depth: 4.25"  
Blade Width: 1"



ISO 9001:2008 • ISO 13485:2003

**FREE TRIAL ON MOST INSTRUMENTS**

Scan to  
Launch Our  
Website



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
Fax +41 41 740 67 71

**INNOMED** The Innomed logo features the company name in a bold, sans-serif font, followed by a stylized lowercase 'i' inside a circle, and a small American flag icon to the right.

**1.800.548.2362**