


# OrthoLucent™ Modified Fukuda-type Retractors

Used to retract the humeral shaft posteriorly, helping to expose the entire glenoid surface, the carbon fiber PEEK composite material is strong, lightweight, completely radiolucent, helps to prevent from marring component surfaces, and can be steam sterilized

PRODUCT NO'S:
1930-R [OrthoLucent™ Narrow] Blade Width: 32 mm Opening: 25 x 40 mm Overall Length: 7.5" (19,1 cm)
1940-R [OrthoLucent™ Wide] Blade Width: 38 mm Opening: 32 x 40 mm Overall Length: 7.5" (19,1 cm)

MADE EXCLUSIVELY  
FOR INNOMED IN  
SWITZERLAND



OrthoLucent™ Wide



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe

Tel. +41 41 740 67 74  
Fax +41 41 740 67 71

1.800.548.2362

# INNOMED