

Offset Femoral Rasp

Designed by Richard Pelliccio

The deep offset design allows the surgeon to line up with canal entry and the tip angled slightly upwards to help prevent femoral protrusion



PRODUCT NO:

4988

Overall Length: 13.67" (34,7 cm)

Handle Length: 5.88" (15 cm)

Rasp Offset: 1.5" (3,8 cm)



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

www.innomed.net

INNOMED 



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362

© 2023 Innomed, Inc.