



Offset Medial Acetabular Retractors with Large Handle

Designed for acetabular exposure during total hip surgery



PRODUCT NO'S:	
C1017 [Left]	
Overall Length: 15.5" (39,4 cm)	
Handle Length: 6" (15,2 cm)	
Blade Width: 32 mm	
C1018 [Right]	
Overall Length: 15.5" (39,4 cm)	
Handle Length: 6" (15,2 cm)	
Blade Width: 32 mm	



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

INNOMED



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362