


Ortho Self-Retaining Retractor with Pin Guides

Designed by Sean Dunn, DPM

Designed to distract a small joint during fusion or osteotomy alignment surgery



- ▶ Parallel pin guides allow pins up to 2 mm
- ▶ Serrated outside blades extend .4" (1 cm) beyond end of pin guides
- ▶ Calibrated ratchet measures the internal diameter of the jaw up to 50 mm

PRODUCT NO:

1842-02

Overall Length: 6.5" (16,5 cm)

Blade Width: 7 mm

Blade Extension (beyond guides): .4" (1 cm)

Blade Thickness: 1.68 mm

Pin Guide Length: 1.25" (3,2 cm)

Pin Guide Internal Diameter: .085" (2,1 mm)



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



INNOMED 



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362