



# OrthoVise™ - Large


U.S. Patent #D398,208

See Reverse Side for  
Small OrthoVise™ Information

- ▶ Made of stainless steel
- ▶ Models equipped with attachment bolts allow a slap hammer to be attached to the end, as well as to either side of the large OrthoVise™ (except Bent Jaw models), for greater adaptability
- ▶ Bent Jaw models are not available with side attachment bolts, but have an end attachment nut to accept a Standard Slap Hammer (#3925 or #3926)
- ▶ A different size slap hammer is used for the large and small sizes of OrthoVise™
- ▶ Slap Hammers are designed with a hammer plate for the additional use of a mallet if desired

## Side Attachment Bolts

Allows a Large OrthoVise™ Slap Hammer (#3950) to be attached to the side of the device. Available on Standard Large and Long Nose Large OrthoVise™ with Attachment Bolts only. (Cannot be added on later.)



## STANDARD LARGE

### PRODUCT NO'S:

	OrthoVise™ Length: 10" (25,4 cm)
3980	with Attachment Bolts (two sides & end) with Large OrthoVise™ Slap Hammer (#3950)
3980-01	with Attachment Bolts (two sides & end) without Slap Hammer
3981	without Attachment Bolts without Slap Hammer with End Attachment Nut that accepts a Standard Slap Hammer (#3925 or #3926)



## LONG NOSE LARGE

### PRODUCT NO'S:

	OrthoVise™ Length: 12" (30,5 cm)
3965	with Attachment Bolts (two sides & end) with Large OrthoVise™ Slap Hammer (#3950)
3965-01	with Attachment Bolts (two sides & end) without Slap Hammer



## LONG NOSE LARGE BENT JAW

### PRODUCT NO'S:

	OrthoVise™ Length: 11.5" (29,2 cm)
3966	with Attachment Nut (end) with Standard Slap Hammer (#3925)
3966-01	without Slap Hammer with Attachment Nut (end) that accepts a Standard Slap Hammer (#3925 or #3926)



## SLAP HAMMERS FOR LARGE ORTHOVISE™ MODELS

### PRODUCT NO'S:

3950	[Slap Hammer for Large OrthoVise] For use with 3965's, 3980's, 3981 Overall Length: 16.5" (41,9 cm)
3925	[Standard Slap Hammer w/16" Rod] For use with: 3966's Overall Length: 16" (40,7 cm)
3926	[Easy Grip Slap Hammer w/16" Rod] For use with: 3966's Overall Length: 16" (40,7 cm)

For Large  
OrthoVise



Standard  
with 16" Rod



Easy Grip Standard  
with 16" Rod



## THREADED ADAPTERS FOR LARGE ORTHOVISE™ MODELS

### PRODUCT NO'S:

3980-02	[Small Adapter] Changes Male End of a Slap Hammer to Female
3980-03	[Threaded Adapting Screw-Large] For use with 3965's, 3966's, 3980's, 3981

Small Adapter



Female/Female  
Adapter Converts  
from Male/Male



Allows a Standard Slap Hammer (#3925 or #3926)  
to be used with any Large OrthoVise™ with Attachment Bolts

Threaded Large



Can be used to append the Large OrthoVise™ with an Attachment Bolt  
for use with a Large OrthoVise™ Slap Hammer (#3950)

MADE EXCLUSIVELY  
FOR INNOMED IN  
GERMANY

[www.innomed.net](http://www.innomed.net)

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
info@innomed-europe.com

1.800.548.2362

# INNOMED



# OrthoVise™ – Small

U.S. Patent #D398,208

See Reverse Side for  
Large OrthoVise™ Information

- ▶ Made of stainless steel
- ▶ Models equipped with attachment bolts allow a slap hammer to be attached to the end of the OrthoVise™ for greater adaptability
- ▶ Slap Hammer is designed with a hammer plate for the additional use of a mallet if desired
- ▶ A different size slap hammer is used for the large and small sizes of OrthoVise™

## STANDARD SMALL

PRODUCT NO'S:	
	OrthoVise™ Length: 8" (20,3 cm)
3985	without Attachment Bolt without Slap Hammer
3985-01	with Attachment Bolt (end) with Small OrthoVise™ Slap Hammer (#3955)
3985-T	with Attachment Bolt (end) without Slap Hammer

## LONG NOSE SMALL


PRODUCT NO'S:	
	OrthoVise™ Length: 9,5" (24,1 cm)
3975	without Attachment Bolt without Slap Hammer
3975-01	with Attachment Bolt (end) with Small OrthoVise™ Slap Hammer (#3955)
3975-T	with Attachment Bolt (end) without Slap Hammer

## SLAP HAMMER FOR SMALL ORTHOVISE™ MODELS

PRODUCT NO:	
3955	[Slap Hammer for Small OrthoVise] For use with: 3975's, 3985's Overall Length: 8.75" (22,2 cm)

## THREADED ADAPTER FOR SMALL ORTHOVISE™ MODELS

PRODUCT NO:	
3985-03	[Threaded Adapting Screw–Small] For use with: 3975's, 3985's



Threaded Small  Can be used to append the small OrthoVise™ with an Attachment Bolt for use with a Small OrthoVise™ Slap Hammer (#3955)

MADE EXCLUSIVELY  
FOR INNOMED IN  
GERMANY



[www.innomed.net](http://www.innomed.net)

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

**INNOMED** 



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
info@innomed-europe.com

1.800.548.2362

© 2023 Innomed, Inc.