

Ott Fat Pad Retractor Set

Designed by Robert Wubben, MD, with modification by David Ott, MD

Designed for soft tissue retraction, the reduced phalange allows for ease of placement in the lateral gutter, and helps avoid contact with the lateral condyle



PRODUCT NO'S:
3219-00 [Set]
Also available individually:
3219-L [Left] Overall Length: 10" (25,4 cm) Prong Length: 20 mm
3219-R [Right] Overall Length: 10" (25,4 cm) Prong Length: 20 mm



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



INNOMED



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362