

# Periosteal Elevator

*Designed for better control*

Designed with a curved end for easier use, and sharper sides for ease of elevating and stripping. The handle is designed for better control.



**PRODUCT NO'S:**

3450 [Curved]  
Overall Length: 7.5" (19,1 cm)  
Handle Length: 4.5" (11,4 cm)  
Blade Size: 16 x 13 mm

3455 [Straight]  
Overall Length: 7.75" (19,7 cm)  
Handle Length: 4.5" (11,4 cm)  
Blade Size: 19 x 14 mm



[www.innomed.net](http://www.innomed.net)

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

**INNOMED** 



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
info@innomed-europe.com

1.800.548.2362