

Pin Puller - Small

Small size allows for use in a small incision to help with removal of a 2 mm or smaller k-wire pin



PRODUCT NO:

3033

Overall Length: 6.5" (16,5 cm)

Jaw Width: 6,2 mm tapering to 3 mm at end

Jaw Height: 11,7 mm



ISO 9001:2008 • ISO 13485:2003

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

INNOMED

1.800.548.2362