



Radiolucent Mini Hohmann Retractors

Designed by Jeffrey Lawton, MD

Radiolucent, lightweight retractors

Carbon fiber material is strong, lightweight, completely radiolucent, can be steam sterilized, and helps to prevent from marring component surfaces.



PRODUCT NO'S:
1591-R [6mm Blade, Bent] Overall Length: 7" Blade Width: 8mm
1592-R [8mm Blade, Deep Bent] Overall Length: 7" Blade Width: 8mm
1593-R [8mm Blade, Bent] Overall Length: 7" Blade Width: 8mm
1594-R [8mm Blade] Overall Length: 6.875" Blade Width: 8mm
1595-R [6mm Blade] Overall Length: 6.875" Blade Width: 6mm
1597-R [16mm Blade] Overall Length: 6.875" Blade Width: 16mm
1596-R [8" Extender] Overall Length: 8"

16mm Blade shown
with Extender



ISO 9001:2008 • ISO 13485:2003

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71



1.800.548.2362