



Radiolucent Richardson-type Soft Tissue Retractors

Designed by Sean B. Kaminsky, MD

Radiolucent, lightweight retractor helps to retract soft tissues for enhanced exposure

Carbon fiber material is strong, lightweight, completely radiolucent, can be steam sterilized, and helps to prevent from marring component surfaces.

Four sizes available.



PRODUCT NO'S:	
3231-23R [23mm] Overall Length: 13" Large Blade: 23mm x 36mm	3231-37R [37mm] Overall Length: 13" Large Blade: 37mm x 52mm
3231-30R [30mm] Overall Length: 13" Large Blade: 30mm x 42mm	3231-44R [44mm] Overall Length: 13" Large Blade: 44mm x 78mm

**PROUDLY
MADE
IN THE
USA**

ISO 9001:2008 • ISO 13485:2003

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

INNOMED 

1.800.548.2362