


# Resnick Allis Bone Clamp

Designed by Charles T. Resnick MD

*A traditional Allis Bone Clamp designed with a longer ratchet which allows for a wider opening to allow a bone to be clamped and locked onto*



**PRODUCT NO:**

1385

Overall Length: 6" (15,2 cm)

Ratcheted Clamp Opens to: 37 mm

Clamp End Width: 4.7 mm



ISO 13485:2016

**FREE TRIAL ON MOST INSTRUMENTS**

[www.innomed.net](http://www.innomed.net)

**INNOMED** 



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
info@innomed-europe.com

**1.800.548.2362**

© 2023 Innomed, Inc.