


Resnick Allis Bone Clamp

Designed by Charles T. Resnick MD

A traditional Allis Bone Clamp designed with a longer ratchet which allows for a wider opening to allow a bone to be clamped and locked onto



PRODUCT NO:

1385

Overall Length: 6" (15,2 cm)

Ratcheted Clamp Opens to: 37 mm

Clamp End Width: 4.7 mm

**PROUDLY
MADE
IN THE
USA**

ISO 9001:2008 • ISO 13485:2003

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

INNOMED 

1.800.548.2362