

Rogozinski Offset Ring Curette

Curette is angled to allow access for easier bone preparation for anterior lumbar interbody fusion

Designed by Chaim Rogozinski, MD and Abe Rogozinski, MD

PRODUCT NO:

5020
Overall Length: 16.25"
Handle Length: 9"
Ring Outside Dia.: 10mm



FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax



INNOMED

1.800.548.2362