

# Romanelli Deep Gelpi Retractor

Designed by Ron Romanelli, MD

*Offers the versatility and ability to be used on minimally invasive total hip replacements*

Can be used to hold the hip capsule out of the way for easy visualization, and to allow reaming of the acetabulum without catching the capsule in the reamer. The ends of the retractor have dull tips to help avoid soft tissue damage. Wider separation occurs at the deep capsule level. Holds the muscle out of the way while retracting the capsule.

**PRODUCT NO:**


4270

Overall Length: 10" (25,4cm)

Depth from Bend: 5.5" (14cm)

Prong Length: 10mm

**PROUDLY  
MADE  
IN THE  
USA**



ISO 9001:2008 • ISO 13485:2003

**FREE TRIAL ON MOST INSTRUMENTS**

Scan to  
Launch Our  
Website



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

**Innomed-Europe** Tel. +41 41 740 67 74  
Fax +41 41 740 67 71

**INNOMED**

**1.800.548.2362**